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Amobarbital Diphenhydramine 3-β-D-glucuronide Quinacrine Amoxicillin 5,5 - Diphenylhydantoin Morphine sulfate Quinidine Ampicillin Doxylamine Nalidixic acid Quinine L-Ascorbic acid Ecgonine hydrochloride Naloxone Ranitidine D-Amphetamine Ecgonine methylester Naltrexone Salicylic acid D,L-Amphetamine sulfate (-)-ψ-Ephedrine Naproxen Secobarbital L-Amphetamine [1R,2S](-) Ephedrine Niacinamide Serotonin Apomorphine L - Epinephrine Nifedipine (5-Hydroxytryptamine) Aspartame Erythromycin Nimesulide Sulfamethazine Atropine β-Estradiol Norcodein Sulindac Benzilic acid Estrone-3-sulfate Norethindrone Sustiva Benzoic acid Ethyl-p-aminobenzoate D-Norpropoxyphene Temazepam Benzoylcgonine Fenoprofen Noscipine Tetracycline Benzphetamine Fursemide DL-Octopamine Tetrahydrocortisone, Bilirubin Gentisic acid Oxalic acid 3- Acetate (±) - Brompheniramine Hemoglobin Oxazepam Tetrahydrocortisone Buspiron Hydralazine Oxolinic acid 3-(β-D glucuronide) Caffeine Hydrochlorothiazide Oxycodone Tetrahydrozoline Cannabidiol Hydrocodone Oxymetazoline Thebaine Cannabinol Hydrocortisone Papaverine Theophyline Chloralhydrate O-Hydroxyhippuric acid Penicillin-G Thiamine Chloramphenicol p-Hydroxyamphetamine Pentazocine Trans-2- Chlordiazepoxide p-Hydroxy- hydrochloride phenylcyclopropylamine Chlorothiazide methamphetamine Pentobarbital Thioridazine (±) - Chlorpheniramine 3-Hydroxytryptamine Perphenazine Tolbutamide Chlorpromazine Imipramine Phencyclidine Trazodone Chlorquine Iproniazid Phenelzine D,L-Tyrosine Cholesterol (±) - Isoproterenol Phenobarbital Triamterene Clomipramine Isosuprine Phentermine Trifluoperazine Clonidine Ketamine Trans-2-phenyl Trimethoprim Cocaethylene Ketoprofen cyclopropylamine Trimipramine Cocaine hydrochloride Labetalol hydrochloride Tryptamine Codeine Levorphanol L-Phenylephrine D,L-Tryptophan Cortisone Loperamide β-Phenylethylamine Tyramine (-) Cotinine Maprotiline Phenylpropanolamine Uric acid Creatinine Meperidine Prednisolone Verapamil Deoxycorticosterone Mephentermine Prednisone Zomepirac.

Methylenedioxy-methamphetamine Concentration (ng/ml)	nper site	Site A		Site B		Site C	
		-	+	-	+	-	+
0	10	10	0	10	0	10	0
250	10	10	0	10	0	10	0
375	10	8	2	9	1	9	1
625	10	1	9	1	9	1	9
750	10	0	10	0	10	0	10

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Acetamidophenol Dextromethorphan Meprobamate Procaine Acetophenetidin Didofenac Methamphetamine Promazine N-Acetylprocainamide Diazepam Methadone Promethazine Acetylsalicylic acid Diflunisal Methoxyphenamine DL-Propranolol Aminopyrine Digoxin Methylphenidate D-Propoxyphene Amitriptyline Dicylamine Morphine- D-Pseudoephedrine

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Methylenedioxy-methamphetamine Concentration (ng/ml)	Percent of Cut-off	n	Visual Result	
			Positive	Negative
0	0	30	30	0
250	-50%	30	30	0
375	-25%	30	25	5
500	Cut-off	30	14	16
625	+25%	30	4	26
750	+50%	30	0	30
1500	3X	30	0	30

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Compound	Concentration (ng/ml)
3,4-Methylenedioxyamphetamine HCl (MDMA)	500
3,4-Methylenedioxyamphetamine HCl (MDA)	3,000
3,4-Methylenedioxyethyl-amphetamine (MDE)	300

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Method	Other MDMA Rapid Test		Total Result
	Positive	Negative	
The MDMA Rapid Test	48	0	48
Dipstick	0	62	62
Total Result	48	62	110
%Agreement With this Rapid Test	>99.9%	>99.9%	>99.9%

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Method	GC/MS		Total Result
	Positive	Negative	
The MDMA Rapid Test	102	1	103
Dipstick	2	145	147
Total Result	104	146	250
%Agreement With GC/MS	98.0%	99.3%	98.8%

1. Winger G. A Handbook of Drug and Alcohol Abuse. Third Edition, Oxford Press. 1992; 146 2. Baselt RC. Disposition of Toxic Drugs and Chemicals in Man. 2nd Ed. Biomedical Publ., Davis, CA. 1982; 488 3. Hawks RL, Chiang CN. Urine Testing for Drugs of Abuse. National Institute for Drug Abuse (NIDA), Research Monograph 73, 1986

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